

Medix Biochemica

Product Manual Cat. No: #2001

Taq 2x PCR Master Mix

Description

The Taq 2x PCR Master Mix contains all the components necessary for PCR, including a Taq DNA polymerase variant and an optimized buffer including ultrapure dNTPs.

Taq 2x PCR Master Mix is a ready to use reaction mix. It contains all components necessary for a successful and reliable PCR or primer extension reaction in all standard PCR cyclers. Only primers and template need to be added.

This mix provides robust PCR performance for a wide range of PCR applications. The pre-ready 2x mix ensures reproducible results, significantly reduces set-up times and the risk of pipetting mistakes.

Applications include standard PCR, real-time-PCR (addition of suitable dye or probe required), primer extension reactions, TA cloning, 3'A-tailing of blunt ends, and screening / high-throughput PCRs.

Kit components

Component	M pack*
Taq 2x PCR Master Mix	1x 1.25 mL

*Other pack sizes, bulk orders and customization are available upon request.

Storage and shipment

Transport with cool packs. The reagents should be stored at -20°C upon arrival. The reagents are stable until the expiration date if stored correctly.

Reaction Master Mix set-up

The recommended master mix set-up for a 50 µL reaction volume is shown in the table below.

Reagent	Volume (µL)	Final concentration
Taq 2x PCR Master Mix	25	1x
∞Forward primer (10 µM)	1	0.2 µM (0.05–1 µM)
∞Reverse primer (10 µM)	1	0.2 µM (0.05–1 µM)
Template / Sample extract	X	<1000 ng* DNA
Nuclease-free water	Up to 50 µL final volume	

Keep all components on ice.

Spin down and mix all solutions carefully before use.

∞Primers should ideally have a GC content of 40–60% typically.

*Suggested template concentration should be about 1 ng – 1000 ng (genomic DNA) or 1 pg – 1 ng (plasmid/viral DNA) per reaction.

Instrument and program set-up

Cycles	Steps	Temperature	Time
1	Initial denaturation	95°C	2 min
25–40	Denaturation	95°C	15 sec
	Annealing*	54–72°C	30 sec
	Extension	72°C	1 min /1000 bp

*Typically, the annealing temperature is about 3–5°C below the calculated melting temperature of the primers used.

Medix Biochemica

Product Manual
Cat. No: #2001