REF 24080

RedStar™ 1kb DNA Ladder

KEY FEATURE

· Suitable for fluorescent-free agarose gel electrophoresis

 $\sqrt{\Sigma}$ 0.5ml

- Contained RedStar[™] DNA/RNA Loading Dye
- Distinct band resolution at 595nm wavelength
- Ideal for 250 bp ~ 10 kb DNA size range

GENERAL USE

- No DNase and RNase detected.
- Load 5 µl per each well of Agarose gel.

For research purpose only. Not for use in diagnostic procedures for clinical purposes.



DESCRIPTION

iNtRON supplies a wide range of products for accurate size and mass estimations (quantitation) of nucleic acid fragments. Nucleic Acid Markers are available for sizing linear, or supercoiled DNA and single-stranded RNA fragments. A variety of these markers are available in the ready-to-use RedStarTM formats.

RedStarTM 1kb DNA Ladder are ideal for determining the size of double-stranded DNA from 250~10,000bp base pairs. The RedStarTM 1kb DNA Ladder consist of 12 linear double-stranded DNA fragments. Several fragments are present at increased intensity to allow easy identification. All fragments are precisely quantified and mixed during the production. For 5 µl loading, all fragments except typical band DNA fragments are 40 ng. The typical band of DNA fragments is 100 ng. These ladders are pre-mixed with loading dye and are ready to use.

All DNA Markers were already stained with RedStar™ DNA/RNA Loading Dye (6x).

STORAGE

- · Store at 4 °C and stable for more than 6 months.
- Store at -20 °C and stable for more than 12 months

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only.

NOTICE BEFORE USE

- · Do not heat before loading.
- Do not use gels containing ethidium bromide (EtBr) or other fluorescent dyes.
- For quantification, adjust the concentration of the sample to equalize, it approximately with the amount of DNA in the nearest band of the ladder.


ELECTROPHORESIS

- · Used for agarose gel without fluorescent material.
- Use 5ul per well ; Usage can be increased or decreased depending on the size of gel and comb.

ORDERING INFORMATIONS

	Product Name	Cat.No.
	RedStar™ DNA/RNA Loading Dye (6x)	21142
	RedStar™ 100 bp DNA Ladder	24079
	RedStar™ 1kb DNA Ladder	24080

DATA EXAMPLES



Size range (bp)	250-10,000	
Conc. (ng/µl)	120	
Typical bands	100ng/5µl	
Other bands	40ng/5µl	
Loading Vol.	5µl	
Band number	12	
Contents	250, 500, 750, 1000,1500, <u>2000,</u> 3000,4000, <u>5000</u> , 6000,8000, 10000	

www.bocascientific.com • (781) 686-1631 • info@bocascientific.com