



Washing Buffer 10x

Wash Buffer for immunoassays

Available products:	Washing Buffer PBS 10x	(article no. 140)
	Washing Buffer PBS 10x without Tween	(article no. 141)
	Washing Buffer TRIS 10x	(article no. 145)
	Washing Buffer TRIS 10x without Tween	(article no. 146)

Storage: 2 – 8 °C or -15 to -30 °C
(tolerates repeated freezing and thawing cycles)

pH-value at 19.0 – 21.0 °C:
(1:10 dilution) 7.2 ± 0.2 (article no. 140, 145, 146)
7.3 ± 0.2 (article no. 141)

Preservative: contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)

Expiry date when stored unopened: see label on the bottle
Use working solution immediately!

For general laboratory use

Instructions for use

Salt crystals may precipitate due to storage at 2 - 8 °C or freezing. Therefore, to prepare the working solution, the *Washing Buffer* must be warmed to room temperature beforehand, which will dissolve any salts that may have precipitated. The buffer should be thoroughly mixed again by shaking immediately before use. The working solution is prepared from the stock solution by diluting 1:10 with distilled or demineralized water and should be used on the same day.

Washing Buffer can be used in automated washers depending on the salt tolerance of the washer. It is essential to observe the specifications of the respective manufacturer!

Suitability of *Washing Buffer* for a specific assay has to be tested by the user.

For further information please visit www.candor-bioscience.com.