Kex2 Protease (Lys/Arg-Arg), recombinant

Product Information Sheet

EP0410/EP0450





SUMMARY

shipped on dry ice; store at -80 °C

For research use only

Product Description and Application

- specific serine endoproteinase
- cleaves amino acid sequence N-Arg-Arg/-C and N-Lys-Arg/-C at the carboxyl end
- suited for protein sequencing and cleavage of fusion proteins with an appropriate recognition sequence

In the series of endoproteinases we also offer the newly developed endoproteinase Lys/Arg-Arg. The specific serine endoproteinase (universal serine-protease) has a MW of 68 kDa. It cleaves at the carboxyl end of the recognition sequences: Arg-Arg/X and Lys-Arg/X and thus provides new possibilities for e.g. the development of fusion protein systems with alternative cleavage sites.

Technical Details

Source: Saccharomyces cerevisiae

Activators: Endoproteinase Lys/Arg-Arg is strongly Ca²⁺ dependent and has a

pH- optimum at pH 7.0

Inhibitors: Ala-Lys-Arg-chloro-methyl ketone

Note: Endoproteinase Lys/Arg-Arg is not inhibited by phenyl-methyl-

sulfonyl-fluoride and tosyl-lysine-chloro-methylketone.

Reaction buffer: 50 mM Tris-HCl or HEPES, 5 mM CaCl₂, pH 7.0 (0.5 mM PMSF; 0.1%

Triton X-100 is not included, but might be added if required by your application). For fusion proteins the reaction conditions have to be

determined empirically.

Substrate solution: Benzyloxycarbonyl-L-tyrosyl-L-lysyl-L-arginin-4-nitroanilid

(Z-L-Tyr-Lys-Arg-pNA \circ TFA salt, Bachem No. L1250). The substrate is dissolved in reaction buffer (OD₃₁₅ = 12). This solution can be

stored at -20 °C for several months. Thaw at room temperature

immediately before use.

Unit definition: 1 unit endoproteinase Lys/Arg-Arg releases 1 µmole 4-nitroaniline per

minute in reaction buffer at pH 7.0 at RT.

Specific activity: 2.2 U/ml

Revised June 2017

Kex2 Protease (Lys/Arg-Arg), recombinant

Product Information Sheet # EP0410/EP0450



Assays:

The substrate solution ($OD_{315} = 12$) was thawed at RT. Assay solution: 100 µl of the substrate solution was mixed with 400 µl reaction buffer and put into a 500 µl cuvette (1 cm optical path length). For this solution the photometer was adjusted to 0.0 and a range from 0 to 0.5 OD_{405} was chosen. The assay solution was mixed with 10 µl of a 1:10 dilution of the endoproteinase Lys/Arg-Arg. The optical density was measured at 405 nm for 5 minutes (at RT). The initial slope was plotted as a straight line and, after 1 minute, the value of the straight line was read (ΔOD_{405} , for 1 minute).

 $\Delta OD_{405} \times 1 \text{ min } \times 100 \times 10$ $= \text{activity in units/min } \times \text{ml (undiluted)}$ $\epsilon \times 2$ (Note: $\epsilon = 9.94$)

References

Achstetter T., Biol. Chem., 259 (1984) 13334-48

Achstetter T., EMBO, 4, No. 1 (1985) 173-177

Dmochowska A., Cell, 50 (1987) 573-584

Fuller R. et al., Arch. of Biochem. and Biophys., 207 No. 2 (1981) 445-454

Fuller R. et al., Arch. of Biochem. and Biophys., 226 (1983) 292-305

Fuller R. et al., Ann. Rev. Physiol., 50 (1988) 345-362

Fuller R., Proc. Acad. Sci. USA, 86 (1989) 1434

Fuller R. et al., Science, 246 (1989) 482

Julius, Cell, 36 (1984) 309

Julius, Cell, 37 (1984) 1075-1089

Mizuno, Biochem. and Biophys. Res. Commun., 144 No. 2 (1987) 807

Wagner J.C., FEBS-Letters, 218 (1987) 31-34

Wagner J.C., FEBS-Letters, 221 No. 2 (1987) 423-426

Revised June 2017

Kex2 Protease (Lys/Arg-Arg), recombinant

Product Information Sheet # EP0410/EP0450



Order Information, Shipping and Storage

Order#	Product	Quantity
EP0410	Kex2 Protease (Lys/Arg-Arg), recombinant	10 U
EP0450	Kex2 Protease (Lys/Arg-Arg), recombinant	50 U
shipped on dry ice; store at -80 °C		

Revised June 2017