

Pro-DeliverIN CRISPR transfection reagent



Protocol

For CRISPR/ CAS9 - Genome Editing using Cas9 protein or Cas9/gRNA RNP complexe.

We Bring The World Of Biotechnology To You

ProDeliverIN CRISPR Reagent User Guide

Package contents	PIC60100: 100µl of ProDeliverIN CRISPR reagent PIC60500: 500µl of ProDeliverIN CRISPR reagent Pro-DeliverIN™ CRISPR is provided with 100µL of R-Phycoerythrin Positive Control.		
Storage Conditions	Store at 4°C upon receipt		
Product Description	ProDeliverIN CRISPR kit is a transfection reagent optimized for recombinant Cas9 protein delivery or Cas9/gRNA RNP complexes. For your gene editing applications, this reagent provides high transfection efficiency with minimal toxicity.		
Important notice	For Research use only. Not for use in diagnostic procedures		

BEFORE YOU BEGIN:

- Protein or Protein/RNP complex solution and Pro-DeliverIN CRISPR should be used at room temperature and gently vortexed prior to use.
- All the complexes must be prepared in medium without serum and supplement.
- It is not recommended to use RPMI during complex preparation, prefer DMEM or PBS.
- For sensitive cells, medium can be replaced with fresh complete culture medium 4 to 6h after transfection.

Tissue Culture Dish	Cell Number per well	Protein or Protein/ RNP complex quan- tity per well*	Total transfection volume per well
96-well 24-well 12-well 6-well 60 mm dish 90-100 mm dish T75 flask	$0.5 - 2.0 \times 10^{4}$ $0.5 - 1.0 \times 10^{5}$ $1.0 - 2.0 \times 10^{5}$ $2.0 - 4.0 \times 10^{5}$ $0.5 - 1.0 \times 10^{6}$ $1.0 - 2.0 \times 10^{6}$ $2.0 - 5.0 \times 10^{6}$	0.4 µg 1 µg 2 µg 5 µg 10 µg 30 µg 35 µg	0.2 mL 0.5 mL 1.0 mL 2.0 mL 4.0 mL 8.0 mL

Table 1: Recommended cell number, protein or protein/RNP complex and transfection volume per well.

PROTOCOL STEPS

The following protocol is given for a single well of a 24-well tissue culture plate containing $^{\sim}1\times10^{5}$ cells/well in 400 μL complete culture serum. If a different culture plate format is used, adjust cell number and reagent amounts according to the table 1.

NOTES:

- -> Pro-DeliverIN CRISPR (PI) should be stored at +4°C. Use 2 μ L of PI per μ g of Cas9 endonuclease or Cas9/gRNA RNP complex.
- -> To find the ideal conditions, Pro-DeliverIN CRISPR must be tested at ratios 1 μ L/ μ g protein, 2 μ L/ μ g protein and 2.5 μ L/ μ g protein.
- -> For the protein quantity, we suggest 0.4 μg per well in 96-well, 1 μg per well in 24-well and 5 μg per well in 6-well.
- -> Depending on the properties of your protein (size, charge...) the **amount** used in the test **can be doubled** (i.e in a 24-well plate, 2 μg of protein instead of 1 μg for 2 / 4 / 5 μL of Pro-DeliverIN CRISPR).

1. Cas9 protein solution

Prepare a suspension of Cas9 endonuclease or RNA complexes at 100 μ g/mL in PBS.

2. Pro-DeliverIN CRISPR solution

Add 2 µL of Pro-DeliverIN CRISPR in a new tube.

3. Complex preparation

Mix 10 μ L of protein or protein/RNP complex solution (1 μ g) with Pro-DeliverIN CRISPR solution. Mix gently by carefully pipetting up and down.

Incubate the mixture for 20 min at room temperature.

4. Transfection

Add the complexes dropwise onto the cells and homogenize by gently rocking the plate side to side to ensure a uniform distribution of the mixture. Incubate the cells under your standard culture conditions for 3 to 48 h.

OPTIONAL

Perform a medium change 2 to 6h after transfection. Withdraw the transfection medium and add fresh growth medium

Additional products for CRISPR Cas9 experiments:

- RmesFect CRISPR for mRNA transfection
- PolyMag CRISPR for Genome editing using expression plasmids
- ViroMag CRISPR to enhance transduction efficiency of CRISPR/ Cas9 viruses

Purchaser Notification

Limited License

The purchase of this product grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in section 1, Kit Contents). This reagent is intended for in-house research only by the buyer. Such use is limited to the purposes described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences. Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of this product. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, please contact us. Buyers may end this License at any time by returning the product and documentation to OZ Biosciences, or by destroying all components. Purchasers are advised to contact us with the notification that the product is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s). This document covers entirely the terms of the research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

This product and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the kit components by following proper research laboratory practices.



www.bocascientific.com (781) 686-1631 info@bocascientific.com