

RmesFect for mRNA transfection ± Lullaby for gRNA transfection

For CRISPR/CAS9 – Genome Editing using mRNA

The following protocol is given for a single well from a 24-well tissue culture plate containing $\sim 1 \times 10^5$ cells/well in 400 μL complete culture serum.

If a different culture plate format is used, adjust cell number and reagent amounts according to the table below:

Tissue Culture Dish	Cell Number per well	Total transfection volume per well
96-well	$0.5 - 2.0 \times 10^4$	0.2 mL
24-well	$0.5 - 1.0 \times 10^5$	0.5 mL
12-well	$1.0 - 2.0 \times 10^5$	1.0 mL
6-well	$2.0 - 4.0 \times 10^5$	2.0 mL
60 mm dish	$0.5 - 1.0 \times 10^6$	4.0 mL
90-100 mm dish	$1.0 - 2.0 \times 10^6$	8.0 mL
T75 flask	$2.0 - 5.0 \times 10^6$	12.0 mL

Table 1: recommended cell number and transfection volume.

Key parameter before beginning the procedure:

- The vector (DNA, gRNA, mRNA, protein, viral particles) and transfection reagents solutions should have an ambient temperature and be gently vortexed prior to use.
- All the complexes must be prepared in medium without serum and supplement.
- It is not recommended to use RPMI during complex preparation, prefer DMEM or PBS.
- For sensitive cells, medium can be replaced with fresh complete culture medium 4 to 6h after transfection.

NOTES: RmesFect (**RM**) and Lullaby (**LL**) should be stored respectively at -20°C and $+4^\circ\text{C}$. Use **3 μL** of RM per μg mRNA and **4 μL** of LL per 50 to 125 nM gRNA.

1. OPTION 1: Cas9 mRNA and gRNA are delivered together

- Cas9 mRNA: dilute **0.5 μg** mRNA into **25 μL** DMEM without any supplement
- Short guide RNA: prepare a solution of gRNA for a final concentration of **50-125 nM** in **25 μL** DMEM without any supplement.
- RmesFect solution: dilute **4 μL** RmesFect into **50 μL** DMEM without any supplement.
- Complexes formation: Mix mRNA suspension with gRNA solution, gently pipette up and down several times and add the mix to RmesFect solution.

Incubate the mixture for **20 min** at room temperature.

2. OPTION 2: Cas9 mRNA and gRNA are delivered separately

- Cas9 mRNA: dilute **0.5 μg** mRNA into **25 μL** DMEM without any supplement
- Short guide RNA: prepare a solution of gRNA for a final concentration of **50-125 nM** in **25 μL** DMEM without any supplement.
- RmesFect solution: dilute **1.5 μL** RmesFect into a new tube containing **25 μL** DMEM without supplement..
- Lullaby solution: dilute **4 μL** Lullaby into a new tube containing **25 μL** DMEM without supplement.
- Complexes formation

Add Cas9 mRNA suspension to RmesFect solution, incubate **20 min** at RT.

Add gRNA solution to Lullaby, incubate **20 min** at RT

3. Transfection

Add the complexes dropwise onto the cells and homogenize by gently rocking the plate side to side to ensure a uniform distribution of the mixture.

Incubate the cells under your standard culture conditions for **24 to 72 h**.

Optional: perform a medium change 2 to 6 H after transfection. Withdraw the transfection medium and add fresh growth medium (optionally without antibiotics)

