

RmesFect CRISPR transfection reagent



Protocol

For RNA (mRNA, gRNA) transfection in CRISPR/ CAS9
Genome Editing

We Bring The World Of Biotechnology To You

RmesFect CRISPR Reagent User Guide

Package contents	RMC70500: 500µl of RmesFect CRISPR reagent		
Storage Conditions	Store at -20°C upon receipt		
Product Description	RmesFect™ CRISPR transfection Reagent based on the TEE- technology is specifically designed for mRNA/gRNA transfection with high efficiency and low toxicity. RmesFect™ is efficient in a large variety of cells.		
Important information	For Research use only. Not for use in diagnostic procedures		

BEFORE YOU BEGIN:

- The RNA solution (mRNA, gRNA) and RmesFect CRISPR should be used at room temperature and be gently vortexed prior to use.
- All the complexes must be prepared in medium without serum and supplement.
- It is not recommended to use RPMI during complex preparation, prefer DMEM or PBS.
- For sensitive cells, medium can be replaced with fresh complete culture medium 4 to 6h after transfection.

Tissue Culture	Cell Number per	RNA (mRNA +	Total transfection volume per well
Dish	well	gRNA) quantity	
96-well 24-well 12-well 6-well 60 mm dish 90-100 mm dish	$0.5 - 2.0 \times 10^{4}$ $0.5 - 1.0 \times 10^{5}$ $1.0 - 2.0 \times 10^{5}$ $2.0 - 4.0 \times 10^{5}$ $0.5 - 1.0 \times 10^{6}$ $1.0 - 2.0 \times 10^{6}$ $2.0 - 5.0 \times 10^{6}$	0.25 μg 0.5 μg 1.0 μg 2.0 μg 4.0 μg 8.0 μg 10.0 μg	0.2 mL 0.5 mL 1.0 mL 2.0 mL 4.0 mL 8.0 mL 12.0 mL

Table 1: Recommended cell number, total RNA quantity and transfection volume per well.

PROTOCOL STEPS

The following protocol is given for a single well of a 24-well tissue culture plate containing $^{\sim}1\times10^{5}$ cells/well in 400 μL complete culture serum. If a different culture plate format is used, adjust cell number and reagent amounts according to the table 1.

NOTES: RmesFect CRISPR (RM) should be stored at -20°C. Use 3 µL of RM per µg total RNA.

1. Cas9 mRNA and gRNA solutions

A. Cas9 mRNA:

Dilute 0.5 μg mRNA into 50 μL DMEM without any supplement

B. Short guide RNA & Cas9 mRNA

Prepare a mix of gRNA and Cas9 mRNA in a 0.25:1 to 1:1 ratios (respectively gRNA to Cas9 mRNA) for a final quantity of 0.5 μ g total in 50 μ L DMEM w/o supplement.

2. RmesFect solution

Dilute 1.5 μ L RmesFect into 50 μ L DMEM without any supplement.

3. Complexes preparation

Mix mRNA suspension or gRNA/mRNA solution with RmesFect CRISPR. Incubate the mixture for 20 min at room temperature.

4. Transfection

Add the complexes dropwise onto the cells and homogenize by gently rocking the plate side to side to ensure a uniform distribution of the mixture. Incubate the cells under your standard culture conditions for 6 to 72h.

OPTIONAL:

Perform a medium change 2 to 6h after transfection. Withdraw the transfection medium and add fresh growth medium

Additional products for CRISPR Cas9 experiments:

- PolyMag CRISPR for Genome editing using expression plasmids
- ProDeliverIN CRISPR for Cas9 protein delivery
- ViroMag CRISPR to enhance transduction efficiency of CRISPR/Cas9 viruses

Purchaser Notification

Limited License

The purchase of this product grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in section 1, Kit Contents). This reagent is intended for in-house research only by the buyer. Such use is limited to the purposes described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences. Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of this product. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, please contact us. Buyers may end this License at any time by returning the product and documentation to OZ Biosciences, or by destroying all components. Purchasers are advised to contact us with the notification that the product is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s). This document covers entirely the terms of the research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

This product and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the kit components by following proper research laboratory practices.



www.bocascientific.com (781) 686-1631 info@bocascientific.com