

# Isohelix **Spin+Collect** TM

Cat. No: SC/ST-100 SC-150

The **Spin+Collect** <sup>™</sup> is a novel sample recovery device to increase sample yields from swabs and other materials.

# **Applications**

Increases sample yields from swabs and other materials. Improves yields of DNA, Nucleotides, Protein, Bacteria etc.

# **Benefits**

Yields increase over 25%.

Recovers sample deeply embedded in the matrix.

Requires a single centrifuge step.

Disposable, sterile and packaged aseptically.

Tested for Human DNA contamination.





The Isohelix **Spin+Collect** firmly grips the swab or other material whilst exposing maximal surface area for the first step of elution or lysis of the sample. A 1.5ml microfuge tube is then fitted as a cap and inverted for centrifugation and final release of any sample trapped in the matrix.

Cat. No:

SC/ST-100 100 x **Spin+Collect** devices, 5 per bag with 100 x 1.5ml matched tubes.

SC-150 150 x **Spin+Collect** devices packed 5 per bag.

Cell Projects Ltd has an exclusive license on patented designs from Guys and St Thomas NHS Foundation to commercialise the **Spin+Collect** and is jointly working with the Foundation to extend its application and use.





## Instructions


- Place Spin+Collect in a standard micro tube rack or dry block heater.
- 2. Place swab (or material plus shaft) firmly into the fins in the cups.
- 3. Snap the swab stick at break point (or cut using an inverted micro tube as a guide. [See next picture]).
- 4. Add extraction solution to a maximum volume of 600ul depending on swab size, matrix and final volume needed.



- 5. Snap a 1.5ml micro tube\* onto the **Spin+Collect**.
- \* Important not all 1.5ml micro tubes will fit correctly!



6. Simply invert, place in a microcentrifuge and spin briefly at full speed to harvest the entire sample from the swab into the micro tube.



7. Discard the **Spin+Collect** and swab and then cap the 1.5ml micro tube containing the sample.



### This has been designed solely for research use only

#### **Other Isohelix Products**

#### Isohelix GeneFix<sup>™</sup> Saliva DNA & RNA collection kits

► GeneFix<sup>™</sup> Saliva DNA & RNA Collection Devices for collecting, storing, stabilising and transporting Saliva Samples

#### Isohelix DNA Buccal Swabs.

> High yields, blood alternative, reproducible, easy to use, different formats for various extraction methodologies.

#### Isohelix DNA / RNA Swab Stabilization

Solutions for physical or chemical stabilization of DNA swabs or chemical to preserve DNA Yield and integrity over extended periods

#### **DNA Isolation and Handling kits**

> DNA Isolation kits, fully optimized to complement the Isohelix range of swabs and GeneFix collectors, using either precipitation or spin column methodologies, both maximising yields and purity.

