Anti-Conjugated Dopamine, monoclonal

Product Information Sheet # 1029GE



SUMMARY

shipped on blue ice; store at -20 °C

For research use only.

Not for use in diagnostic or therapeutic procedures.

Description

Monoclonal antibody was obtained after BALB/c mouse immunization with the conjugates: Dopamine-Glutaraldehyde-Carrier proteins and hybridization of spleen cells with the myeloma cell line SP2/O/Ag14. Ascite production was performed in BALB/c mice.

Subclass: IgG 1, Kappa

Purification

The ascitic fluid was purified by ammonium sulfate precipitation and gel filtration.

Contents: The antibody is supplied in a unit size of 50 µl in PBS with 50% glycerol.

Storage and Handling

When this product is stored undiluted at -20 °C or lower, it is stable for at least 2 years. Multiple freezing and thawing is only possible up to 5 times. The antibodies can be prediluted 10x in PBS containing 0.1% merthiolate or a mixture of PBS/glycerol (1:1, vol./vol.). At 4 °C these solutions are stable for 2 months.

Specificity

Using a conjugate Dopamine-Glutaraldehyde-Protein, antibody specificity was determined in an ELISA test by competition experiments with the following compounds:

Compound	Cross-reactivity ratio (a)	
Dopamine-G-BSA	1	
L-DOPA-G-BSA	1:	10,000
Tyrosine-G-BSA	1:	36,000
Tyramine-G-BSA	1:	> 100,000
Noradrenaline-G-BSA	1:	> 100,000
Octopamine-G-BSA	1:	> 100,000
Adrenaline-G-BSA	1:	> 100,000
L-DOPA=G=BSA (b)	1:	> 100,000
Dopamine	1:	> 100,000

(a): Dopamine-G-BSA concentration/unconjugated or conjugated catecholamine concentration at half displacement;

(b): Non-reduced conjugate;

G = Glutaraldehyde, BSA = Bovine Serum Albumine.

MoBiTec GmbH, Germany Phone: +49 551 70722 0 Fax: +49 551 70722 22 E-Mail: info@mobitec.com www.mobitec.com

Anti-Conjugated Dopamine, monoclonal

Product Information Sheet # 1029GE



Recommended dilution

The antibody was tested using the free floating PAP technique on rat dopaminergic areas. The anticonjugated Dopamine antibody gave a good staining between a 1:5,000-1:20,000 dilution in these areas.

References

CHAGNAUD J.L., MONS N., TUFFET S., GRANDIER-VAZEILLES X. and GEFFARD M. Monoclonal antibodies against glutaraldehyde-conjugated dopamine. J. Neurochem., (1987) 49, 487-494.

DECAVEL C., LESCAUDRON L., MONS N. and CALAS A. First visualization of dopaminergic neurons with a monoclonal antibody to dopamine: a light and electron microscopic study. J. Histochem. Cytochem., (1987) 35, 1245-1252.

SCHÜRMANN F.W., ELEKES K. and GEFFARD M. Dopamine-like immunoreactivity in the bee brain. Cell & Tissue Res., (1989), <u>256</u>, 399-410.

GOLDMAN-RAKIC P.S., LERANTH C., WILLIAMS S.M., MONS N. and GEFFARD M. Dopamine synaptic complex with pyramidal neurons in primate cerebral cortex. Proc. Nat. Acad. Sci. USA, (1989), 86, 9015-9019.

Order Information, Shipping and Storage

Order#	Product	Quantity
1029GE	Anti-Conjugated Dopamine, monoclonal	50 μΙ
shipped on blue ice; store at -20 °C		

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders Technical Service – Product information

phone: +49 (0)551 707 22 0 phone: +49 (0)551 707 22 70 fax: +49 (0)551 707 22 27 e-mail: order@mobitec.com e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com

MoBiTec GmbH, Germany Phone: +49 551 70722 0 Fax: +49 551 70722 22 E-Mail: info@mobitec.com www.mobitec.com