

Product Information Sheet # 1076GE

SUMMARY

shipped on blue ice; store at -20 °C

For research use only.

Not for use in diagnostic or therapeutic procedures.

Description

Polyclonal antisera were raised in rabbits after immunization with the bacterial antigen embedded in a 2% agarose gel.

Contents: The antibody is supplied in a unit size of 50 µl in PBS with 50% glycerol.

Storage and Handling

When this product is stored undiluted at -20 °C or lower, it is stable for at least 2 years. Multiple freezing and thawing is only possible up to 5 times. The antibodies can be prediluted 10x in PBS containing 0.1% merthiolate or a mixture of PBS/glycerol (1:1, vol./vol.). At 4 °C these solutions are stable for 2 months.

Specificity

Antibody specificity was determined in an ELISA test by competition experiments with the following compounds:

Concentration

Co	m	р	0	und	
-					

Bacterial antigen from the following strains:	
Hafnia alvei	≥ 10 ⁻² g/l
Morganella morganii	≥ 10 ⁻² g/l
Pseudomonas putida	≥ 10 ⁻² g/l
Proteus mirabilis	≥ 10 ⁻² g/l
Citrobacter koseri	≥ 10 ⁻² g/l
Klebsiella pneumoniae	≥ 10 ⁻² g/l

Affinity for Morganella morganii preparation: 0.16 x 10⁻³ g/l

Recommended dilution

The anti-Pseudomonas aeruginosa antibodies can be diluted 1:5,000 to 1:10,000.

Product Information Sheet # 1076GE



Order Information, Shipping and Storage

Order#	Product	Quantity
1076GE	Polyclonal anti-Pseudomonas aeruginosa antibodies	50 µl
shipped on I	olue ice; store at -20 °C	

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service – General inquiries & orders Technical Service – Product information

phone: +49 (0)551 707 22 0	phone: +49 (0)551 707 22 70
fax: +49 (0)551 707 22 22	fax: +49 (0)551 707 22 77
e-mail: order@mobitec.com	e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com