

## Genomic DNA isolation procedure for amniotic fluid samples

1. Take at least 2 ml of amniotic fluid fresh or 50000-150000 cells of cultured amniotic fluid, previously resuspended, and put it into a provided 2 ml tube.

**Note:** Yield depends on gestation week of the fetal material from which the sample was sourced. The kit has been tested to yield  $>1\mu$ gr using starting materials of approximately 2 ml of amniotic fluid obtained from fetal amniocentesis at 15 weeks of gestation or more. However, using larger volume when possible is recommended.

- 2. Centrifuge at 11200 g (approximately 12000 rpm) for 5 minutes at room temperature.
- 3. Discard the supernatant with care.

Note: Repeat steps 1 to 3 if sample volume is larger than 2 ml, until all sample has been centrifuged.

- Add 300 µl of Solution A and vortex gently.
  Note: Mix gently before addition if Solution A was stored at 4°C
- Incubate at 55°C for 30 minutes. Use of a horizontal shaker (at 100-150 rpm) is optional but preferable. <u>Optional RNase treatment</u>: heat sample at 70°C for 5 minutes. Add RNase (not provided) 100 µg/ml (final concentration) and incubate at 37°C for 5 minutes. Finally, continue to step 6.
- 6. Cool the samples (from step 5) by incubating approximately 5 minutes on ice or refrigerator.
- 7. Add 100 µl of Solution B and vortex gently.
- 8. Centrifuge at 17000 g (approx. 13500 rpm) for 10 min., at room temperature
- 9. Carefully pipet the supernatant into a provided 1.5 ml tube, discarding the remaining pellet.
- 10. Add 40  $\mu l$  of Solution C and 400  $\mu l$  of Solution D.
- 11. Shake slightly by inverting the tube several times until a homogenous solution is observed.
- 12. Incubate the tube at room temperature for 10 minutes in a vertical position.
- Centrifuge at 11200 g (12000rpm) for 5min at room temperature. Note: Usually, a little pellet forms.
- 14. Discard the supernatant with care.
- 15. Add 500 µl of Solution E.
- Centrifuge at 11200 g (12000 rpm) for 5 minutes at room temperature.
  Optional: An additional wash with 70% ethanol can be done to avoid salt excess in the final sample.
- 17. Discard the supernatant with care (\*) and place the tube (containing the pellet) open for 5-10 minutes to dry the pellet.

(\*)Pellet is easily removable from the bottom of the tube, so supernatant must be discarded very carefully to avoid sample loss.

- Add 30-50 µl of Solution F and pipet up and down carefully to resuspend the genomic DNA. Note: Nuclease-free water can be used but is not recommended for long time storage.
- 19. Optional, incubate the tube at 37°C for 30 minutes to help DNA solubilisation.
- 20. Use immediately or store at 4°C (if used during the next 48 hours) or at -20°C (for longer storage).



## Genomic DNA isolation procedure for chorionic villus samples

- 1. Take a small amount of chorionic villi (1-5 mg) or 50000-150000 cells of cultured chorionic villi and place it into a provided 2 ml tube.
- Add 300 µl of Solution A and vortex gently. If you start with cell culture first centrifuge at 11200 g (approximately 12000 rpm) for 5min at room temperature and discard the supernatant with care. Then add 300 µl of Solution A.

Note: Mix gently before addition if Solution A was stored at 4°C.

- Incubate at 55°C for 60 minutes. Use of a horizontal shaker (at 100-150 rpm) is optional but preferable <u>Optional RNase treatment</u>: heat sample at 70°C for 5 minutes. Add RNase (not provided) 100 µg/ml (final concentration) and incubate at 37°C for 5 minutes. Finally, continue to step4.
- 4. Cool the samples (from step 3) by incubating approximately 5 minutes on ice or refrigerator.
- 5. Add 100 µl of Solution B and vortex gently.
- 6. Centrifuge at 17000 g (approx. 13500 rpm) for 10 min., at room temperature.
- 7. Carefully pipet the supernatant into a provided 1.5 ml tube, discarding the remaining pellet.
- 8. Add 40 µl of Solution C and 400 µl of Solution D.
- 9. Shake slightly by inverting the tube several times until a homogenous solution is observed.
- 10. Incubate the tube at room temperature for 10min in a vertical position.
- Centrifuge at 11200 g (12000 rpm) for 5 minutes, at room temperature. Note: Usually a little pellet forms.
- **12.** Discard the supernatant with care.
- 13. Add 500 µl of Solution E.
- **14.** Centrifuge at 11200 g (12000 rpm) for 5 minutes at room temperature.

Optional: An additional wash with 70% ethanol can be done to avoid salt excess in the final sample.

**15.** Discard the supernatant with care(\*) and place the tube (containing the pellet) open for 5-10 minutes to dry the pellet.

(\*)Pellet is easily removable from the bottom of the tube, so supernatant must be discarded very carefully to avoid sample loss.

- Add 30-50 µl of Solution F and pipet up and down carefully to resuspend the genomic DNA.
  Note: Nuclease-free water can be used but is not recommended for long time storage.
- 17. Optional, incubate the tube at 37°C for 30 minutes to help DNA solubilisation.
- 18. Use immediately or store at 4°C (if used during the next 48 hours) or at -20°C (for longer storage).