

Viral Applications

# LentiBlast Premium

Transduction reagent Enhance lentiviral infection and tranduction efficiency







# LentiBlast Premium Quick Protocol

Use the following protocol to find the ideal conditions for LentiBlast Premium in 24-well plate. If the lentiviral transduction/infection conditions are unknown, we recommend starting with a MOI of 2 using a lentiviral vector encoding for a fluorescent protein.

**NOTE:** We suggest to use 0.5 to 10 µL of LentiBlast Premium per conditions. For hard-to-transduce cells, you can add a centrifugation step to this protocol.

## **1.** VIRUS PREPARATION



Dilute virus into culture medium sufficient for 4 samples (50 µL each).

MOI 2 is recommended in case of unknown lentiviral transduction conditions.

## **2.** DISPATCH EQUAL VOLUME OF VIRAL SUSPENSION INTO 4TUBES



## **3.** ADD LENTIBLAST TO EACH TUBE

	Tube 1	Tube 2	Tube 3	Tube 4
LentiBlast	-	0.5 μL	5 μL	10 µL

## 4. MIX VIALS BY INVERTING

Do not vortex or centrifuge

- **5.** ADD VIRUS +/- LENTIBLAST Incubate the cells 24 h under standard culture conditions
- 6. Optional: 24 h MEDIUM EXCHANGE Remove medium from the cells add pre-warmed culture medium

## **7.** INCUBATE CELLS 24 TO 96 h.

Incubate the cells under standard culture conditions We recommend performing assay from 24 to 96 h.







## **IMPORTANT NOTES – Before you begin**

- $\checkmark$  Allow reagents to reach RT and gently vortex them before utilization.
- ✓ Effects of lentiviral transduction are generally observed after 48 to 96 h.
- ✓ For hard-to-transduce cells, you can add a centrifugation step to this protocol.
- ✓ Dilute your reagent with deionized water for doses less than 1µL.
- $\checkmark$  If transduction conditions are unknown, we recommend beginning with a MOI of 2.
- ✓ Do not use LentiBlast Premium with another viral enhancer or adjuvant.
- ✓ In case of low efficiency using the standard protocol, transduction efficiency can be raised:
  - a. **by centrifugation:** However, this procedure is not recommended because in some cases it can lower the viability. Please refer to Centrifugation Protocol below for more details.
  - b. by optimization: Variations in MOI, volumes of LentiBlast Premium can lead to higher transduction efficiency.

Package content	LBPX500: 500 µL of LentiBlast Premium Reagent LBPX1500: 1500 µL of LentiBlast Premium Reagent		
Shipping conditions	Room Temperature		
Storage conditions	Store the LentiBlast Premium Reagent at -20°C upon reception		
Shelf life	1 year from the date of purchase when properly stored and handled		
Product description	LentiBlast Premium is ideal to enhance lentiviral infection and transduction in any type of cells: adherent or in suspension, primary or cell lines. Highly recommended for Stem cells		
Important notice	For research use only. Not for use in diagnostic procedures		

#### 1. Cells Preparation

Cell culture prior to transduction: the day before transduction prepare the cells according to the table below. Cells should be 20-50 % confluent at the time of transduction (see the suggested cell number in the Table 1).

Tissue Culture Dish	Cell Number		
96 wells	3 – 8 x 1.10 <sup>3</sup>		
24 wells	2 – 4 x 1.10 <sup>4</sup>		
6 wells	1 – 2 x 1.10 <sup>5</sup>		

Table 1: Suggested cell number for lentiviral transduction (per well)

### **IMPORTANT NOTE**

For hard-to-transduce or non-permissive cells, prepare the cells the day of transduction and then refer to Centrifugation Protocol.

#### 2. Standard Protocol

Use the quick protocol above to find the ideal conditions for LentiBlast Premium in 24-well plate. If the lentiviral transduction/infection conditions are unknown, we recommend starting with a MOI of 2 using a lentiviral vector encoding for a fluorescent protein.

**NOTE:** We suggest using 0.5, 5  $\mu$ L and 10  $\mu$ L of LentiBlast Premium per conditions.

#### 3. Centrifugation Protocol

For hard-to-transduce cells, it is recommended to add a centrifugation step to the standard protocol. Cells are prepared the day of transduction, counted, pelleted and suspended in Lentivirus/LentiBlast Premium mixes.

NOTE: Centrifugation may influence cell viability.

- 1) Detach cells and seed them into 5 wells. Refer to Table 1 for suggested cell density
- 2) Follow steps 1 to 4 of the quick protocol and add lentivirus/LentiBlast Premium mixes to cells
- 3) Centrifuge the plate 900 rpm for 90 min
- 4) Incubate cells overnight and proceed to steps 7 and 8 of the standard protocol

#### 1. Optimizing LentiBlast Premium volumes

To find the ideal transduction conditions using LentiBlast Premium, we recommend optimizing volumes of LentiBlast Premium with a fixed MOI (refer to Table 2).

Tissue culture dish format	1:1000	1:500	1:250	1:100	1:50	1:25
96-well plate	0.15 µL	0.3 µL	0.6 µL	1.5 µL	3 µL	6 µL
24-well plate	0.5 µL	1 µL	2 µL	5 µL	10 µL	20 µL
12-well plate	1μL	2 µL	4 µL	10 µL	20 µL	40 µL
6-well plate	2 µL	4 µL	8 µL	20 µL	40 µL	80 µL

Table 2: Recommended dilutions and volumes of LentiBlast Premium depending on tissue culture dish format

#### 2. Finely tuning transduction parameters

We have observed that replacing culture medium at the time of transduction with complete medium containing viral particles/LentiBlast Premium could improve the overall efficiency.

- 1) Seed the cells as previously described
- 2) Prepare a viral suspension in complete culture medium; volume should correspond to cell culture volume (refer to Table 1 for recommended volumes)
- 3) Add LentiBlast Premium to the viral suspension
- 4) Mix vial by inverting
- 5) Replace cell culture medium with viral suspension/LentiBlast Premium solution
- 6) Optionally perform a medium change after 24H, and incubate 24H to 96H