


ViroMag CRISPR transduction reagent



Protocol

For CRISPR/CAS9 viral vectors

PolyMag CRISPR Reagent User Guide

Package contents	VMC50200: 200µl of ViroMag CRISPR reagent VMC51000: 1mL of ViroMag CRISPR reagent KVC50100 : 100µL of ViroMag CRISPR reagent + 1 Super Magnetic plate
Storage Conditions	Store at 4°C upon receipt
Product Description	ViroMag CRISPR reagent is the only magnetic viral transduction enhancer for CRISPR/Cas9 viruses (adenovirus, lentivirus, retrovirus...) for gene editing applications.
Important guidelines	This reagent has to be used with a magnetic plate For Research use only. Not for use in diagnostic procedures

BEFORE YOU BEGIN:

- The viral vector solution and ViroMag CRISPR Should be used at room temperature and be gently vortexed prior to use.
- All the complexes must be prepared in medium without serum and supplement.
- It is not recommended to use RPMI during complex preparation, prefer DMEM or PBS.
- For sensitive cells, medium can be replaced with fresh complete culture medium 4 to 6h after transfection.

Tissue Culture Dish	Cell Number per well	ViroMag CRISPR volume per well	Total transduction volume per well
96-well	0.5 – 2.0 × 10 ⁴	1.5 µL	0.2 mL
24-well	0.5 – 1.0 × 10 ⁵	6 µL	0.5 mL
12-well	1.0 – 2.0 × 10 ⁵	12 µL	1.0 mL
6-well	2.0 – 4.0 × 10 ⁵	30 µL	2.0 mL
60 mm dish	0.5 – 1.0 × 10 ⁶	60 µL	4.0 mL
90-100 mm dish	1.0 – 2.0 × 10 ⁶	150 µL	8.0 mL
T75 flask	2.0 – 5.0 × 10 ⁶	150 µL	12.0 mL

Table 1: Recommended cell number, ViroMag CRISPR volumes and transduction volume per well.

PROTOCOL STEPS

The following protocol is given for a single well from a 24-well tissue culture plate containing $\sim 1 \times 10^5$ cells/well in 400 μL complete culture serum. If a different culture plate format is used, adjust cell number and reagent amounts according to the table 1.

NOTES: ViroMag CRISPR should be stored at $+4^\circ\text{C}$. Use 6 μL of ViroMag CRISPR per well. Viral particles will be diluted in 100 μL of medium without any supplement (DMEM) and added onto ViroMag CRISPR magnetic nanoparticles.

1. Viral particle solution

The total MOI may vary depending on viral constructions and cell type: MOI of 3 for HEK-293 cells, and MOI of 60 for iPS. We recommend using a MOI of 10 to begin.

Prepare a viral suspension* for a final MOI=10 in 100 μL of DMEM without any supplement.

2. ViroMag CRISPR solution

Add 6 μL of ViroMag CRISPR in a new tube.

3. Complexes preparation

- A. Add the viral solution onto ViroMag CRISPR nanobeads, mix gently by carefully pipetting up and down.
- B. Incubate the mixture for 20 min. at room temperature.

4. Infection

- A. Add the magnetic complexes dropwise onto the cells and homogenize by gently rocking the plate side to side to ensure a uniform distribution of the mixture.
- B. Place the cells on the magnetic device and incubate at 37°C for 30 min.
- C. After 30 min of incubation, remove the magnetic plate.
- D. Incubate the cells under your standard culture conditions for 24 to 72h.

OPTIONAL:

Perform a medium change. Keep the magnetic plate beneath the cell culture dish, withdraw the transfection medium and add fresh growth medium. Then, remove the magnetic plate.

* Single virus containing both Cas9 & gRNA or multiple viral particles coding for Cas9 and gRNA.

Additional products for CRISPR Cas9 experiments:

- PolyMag CRISPR for Genome editing using expression plasmids
- RmesFect CRISPR for mRNA transfection
- ProDeliverIN CRISPR for Cas9 protein delivery

Purchaser Notification

Limited License

The purchase of this product grants the purchaser a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in section 1, Kit Contents). This reagent is intended for in-house research only by the buyer. Such use is limited to the purposes described in the product manual. In addition, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences. Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of this product. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, please contact us. Buyers may end this License at any time by returning the product and documentation to OZ Biosciences, or by destroying all components. Purchasers are advised to contact us with the notification that the product is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s). This document covers entirely the terms of the research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

This product and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the kit components by following proper research laboratory practices.

We Bring The World Of Biotechnology To You



www.bocascientific.com

(781) 686-1631

info@bocascientific.com