

Product Datasheet

qGreen Mix Separate ROX #2301

Product Name	qGreen Mix Separate ROX
Catalog Number	#2301
Description	qGreen Mix is a universal intercalating dye mix for robust, sensitive, and fast qPCR. qGreen Mix uses state-of-the-art technologies with an antibody-regulated hot-start Taq polymerase and intercalating dye for real-time PCR amplification. The optimized buffer chemistry and PCR enhancers and stabilizers enable rapid and sensitive qPCR.
Tested Applications	End-Point, Real-Time
Brand	Medix Biochemica
Storage	-20°C
Content	100 rxn: 1 x 1 ml qGreen Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive 500 rxn: 5 x 1 ml qGreen Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive