

Product Datasheet

qProbe Mix Separate ROX #2401

Product Name	qProbe Mix Separate ROX
Catalog Number	#2401
Description	qProbe Mix is a universal one-step probe mix for robust, sensitive, and fast qPCR. The mix uses state-of-the-art technologies with an antibody-regulated hot-start Taq DNA polymerase for real-time PCR amplification of single or multiplex DNA targets. qProbe Mix is compatible with several probes such as TaqMan® and Scorpions®.
Tested Applications	End-Point, Real-Time
Brand	Medix Biochemica
Storage	-20°C
Content	100 rxn: 1 x 1 ml qProbe Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive 500 rxn: 5 x 1 ml qProbe Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive