

Product Datasheet

qPCR Lyo-Ready Mix #4401

Product Name	qPCR Lyo-Ready Mix
Catalog Number	#4401
Description	qPCR Lyo-Ready Mix is a universal probe mix that allows robust, sensitive, and fast qPCR for the detection of single or multiplex DNA targets. qPCR Lyo-Ready Mix contains optimized excipients, buffer, PCR enhancers and an antibody-regulated hot-start Taq DNA polymerase and is ready to be lyophilized to produce stable reagents at room temperature. Upon addition of target specific primers/probes to the mastermix, the mixture can be lyophilized directly, without the need to add additional excipients. qPCR Lyo-Ready Mix is compatible with several probes such as TaqMan® and Scorpions®.
Tested Applications	End-Point, Real-Time, Lyophilized
Brand	Medix Biochemica
Storage	-20°C
Content	500 rxn: 4 x 1.25 ml qPCR Lyo-Ready Mix (2x) 10 000 rxn: 100 ml qPCR Lyo-Ready Mix (2x)