

Product Datasheet

qPCR Probe 2 x LyoCake Master Mix (freeze-dried) #9801

Product Name	qPCR Probe 2 x LyoCake Master Mix (freeze-dried)
Catalog Number	#9801
Description	qPCR Probe 2x Lyocake Master Mix is a ready-to-use reaction mix for sensitive and reliable probe-based qPCR in all standard real-time PCR cyclers. It includes an engineered, aptamer-based fast-start DNA polymerase, optimized reaction buffer and ultrapure dNTPs. Freeze-dried qPCR Probe 2x Lyocake Master Mix can be stored at room temperature. For more information, please check : qPCR Probe 2 x LyoCake Master Mix (freeze-dried)
Tested Applications	End-Point, Real-Time, Lyophilized
Brand	myPOLS Biotec
Storage	Ambient / -20°C
References	Aebischer-Gumy, C., Moretti, P., Ollier, R., Ries Fecourt, C., Rousseau, F., & Bertschinger, M. (2020, January). SPLICELECT™: an adaptable cell surface display technology based on alternative splicing allowing the qualitative and quantitative prediction of secreted product at a single-cell level. In MABs (Vol. 12, No. 1, p. 1709333). Taylor & Francis.
Content	S pack: 4 x LyoCakes, 1 x 1.25 ml Rehydration buffer L pack: 20 x LyoCakes, 5 x 1.25 ml Rehydration buffer