Bacillus megaterium protoplasts, strain WH320

MO BI TEC MOLECULAR BIOTECHNOLOGY

Product Information Sheet # BMEG02

SUMMARY shipped on dry ice; store at -80 °C For research use only!



Product Description and Application

• *Bacillus megaterium* protoplasts, strain WH320. Host strain for *Bacillus megaterium* gene expression systems

For handling carefully read the handbook supplied with the gene expression system. Care must be taken not to interrupt the cooling chain. Always store frozen protoplasts at -80 °C or lower.

Genotype

 n/a – WH320 is a chemical mutant of DSM319, lacZ⁻ (Rygus T, Hillen W.: Inducible high-level expression of heterologous genes in *Bacillus megaterium* using the regulatory elements of the xylose-utilization operon, Appl. Microbiol. Biotechnol. (1991), 35, 594-599)

Contents

• 5x500 µl ready-to-use protoplasts

Order Information, Shipping and Storage

Order#	Product	Quantity
BMEG02	Bacillus megaterium protoplasts, strain WH320	5x500 µl
shipped on dr	y ice; store at -80 °C	