

Bacillus megaterium protoplasts, strain YYBm1

Product Information Sheet
BMEG04



SUMMARY

shipped on dry ice; store at -80 °C

For research use only!



Product Description and Application

- *Bacillus megaterium* protoplasts, strain YYBm1. Host strain for *Bacillus megaterium* gene expression systems. Ready for transformation.
- YYBm1 is unable to utilize xylose as carbon source.

For handling carefully read the handbook supplied with the gene expression system. Care must be taken not to interrupt the cold chain. Always store frozen protoplasts at -80 °C or lower.

Genotype

- Defined deletion mutant of wild type strain DSM319: $\Delta nprM$, $\Delta xyIA$

Contents

- 5 x 500 μ l protoplasts of *Bacillus megaterium* YYBm1 supplemented with 11% glycerol.

Order Information, Shipping and Storage

Order#	Product	Quantity
BMEG04	<i>Bacillus megaterium</i> protoplasts, strain YYBm1	5 x 500 μ l
shipped on dry ice; store at -80 °C		