Bacillus megaterium protoplasts, strain MS941

Product Information Sheet # BMEG50



SUMMARY

shipped on dry ice; store at -80 °C

For research use only!



Product Description and Application

• Bacillus megaterium protoplasts, strain MS941. Host strain for Bacillus megaterium gene expression systems

For handling carefully read the handbook supplied with the gene expression system. Care must be taken not to interrupt the cooling chain. Always store frozen protoplasts at -80 °C or lower.

Genotype

• Defined deletion mutant of wildtype strain DSM319, ΔnprM

Contents

Revised February 2017

• 5x500 µl ready-to-use protoplasts

Order Information, Shipping and Storage

Order#	Product	Quantity
BMEG50	Bacillus megaterium protoplasts, strain MS941	5x500 µl
shipped on dry ice; store at -80 °C		