

***Lactococcus lactis* expression strain MG1363**

Product Information Sheet
VS-ELS01363



SUMMARY

shipped on dry ice; store at -80 °C

For research use only!

We Bring The World Of Biotechnology To You



Product Description and Application

- *Lactococcus lactis* subsp. *cremoris* expression strain MG1363. Host strain for *Lactococcus* gene expression systems.

For handling carefully read the handbook pertaining to the gene expression system.
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

Genotype

- Plasmid-free derivative of SH4109 (Gasson M. J.: Plasmid complements of *Streptococcus lactis* NCDO 712 and other lactic streptococci after protoplast-induced curing, Journal of Bacteriology (1983), Vol. 154, No. 1, 1 - 9)

Contents

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

Order Information, Shipping and Storage

Order#	Product	Quantity
VS-ELS01363	<i>Lactococcus lactis</i> expression strain MG1363	1 ml
shipped on dry ice; store at -80 °C		