

# ***Lactococcus lactis* strain NZ3910, food grade**

Product Information Sheet  
# VS-ELS03910-01



## **SUMMARY**

shipped on dry ice; store at -80 °C

For research use only!

*We Bring The World Of Biotechnology To You*



## **Product Description and Application**

- *Lactococcus lactis* strain NZ3910, food grade. Host strain for NICE® *Lactococcus* gene expression systems.

For handling carefully read the handbook pertaining to the gene expression system.  
Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

## **Genotype**

- *lacF*<sup>-</sup>, *nisRnisk* (in neutral locus)

## **Contents**

- 1 ml cell suspension in growth medium supplemented with 10 - 30% glycerol.

## **Order Information, Shipping and Storage**

Order#	Product	Quantity
VS-ELS03910-01	NICE® <i>Lactococcus lactis</i> expression Strain NZ3900, food grade	1 ml
shipped on dry ice; store at -80 °C		