Escherichia coli strain MC1061

Product Information Sheet # VS-ELS10610-01



SUMMARY

shipped on dry ice; store at -80 °C

For research use only!



Product Description and Application

• Escherichia coli strain MC1061. Cloning strain for diverse Gram-positive bacteria, e.g., Lactococcus lactis or Bacillus subtilis.

For handling carefully read the handbook pertaining to the gene expression system. Care must be taken not to interrupt the cooling chain. Always store frozen stock at -80 °C or lower.

Genotype

• araD139, Δ(ara, leu)7697, ΔlacX74, galU-, galK-, hsr-, hsm+, strA (Casadaban MJ and Cohen SN, 1980; Analysis of gene control signals by DNA fusion and cloning in Escherichia coli.)

Contents

• 1 ml cell suspension in growth medium supplemented with 10-30% glycerol.

Order Information, Shipping and Storage

| Order# | Product | Amount |
|-------------------------------------|--------------------------------|--------|
| VS-ELS10610-01 | Escherichia coli strain MC1061 | 1 ml |
| shipped on dry ice; store at -80 °C | | |

www.bocascientific.com • (781) 686-1631 • info@bocascientific.com