

# 10 ml columns (90 µm filter pore size)

Product Information Sheet  
# S10149



We Bring The World Of Biotechnology To You  
**BOCA**  
**SCIENTIFIC**

## SUMMARY


shipped at RT; store at RT

For research use only

Easy-To-Handle Lab Columns with Luer-Lock Connections for Syringes and Tubing.

## Product Description and Application


Set of empty columns for e.g. affinity and gel chromatography, ion exchange columns, gel filtration, etc.



- All columns:**
- Luer-lock connectors at both sides for connecting e.g. syringes or stopcocks
  - compatible with laboratory standard
  - can be closed with additional closing cap
  - smallest and largest volumes can be treated easily
  - come with lower filters with 10 µm, 35 µm or 90 µm pore size inserted
  - upper filter (used on top of the matrix) available separately
  - 2 one-way stopcocks in each package of 20 columns
  - Luer adapters available
  - very competitive price!
- 2.5 ml column:**
- can be closed tightly with Luer-lock closing cap
  - tight up to 4 atmospheres pressure
  - can be centrifuged and autoclaved (with filters at 110 °C)

# 10 ml columns (90 µm filter pore size)

Product Information Sheet  
# S10149



## Order Information, Shipping and Storage

### Laboratory Column

Order#	Product	Quantity
S10121	2.5 ml columns (10 µm filter pore size)	20
S1012	2.5 ml columns (35 µm filter pore size)	20
S10129	2.5 ml columns (90 µm filter pore size)	20
S10131	5.0 ml columns (10 µm filter pore size)	20
S1013	5.0 ml columns (35 µm filter pore size)	20
S10139	5.0 ml columns (90 µm filter pore size)	20
S10141	10 ml columns (10 µm filter pore size)	20
S1014	10 ml columns (35 µm filter pore size)	20
S10149	10 ml columns (90 µm filter pore size)	20
S1210	Upper filter 10 µm for 2.5 ml columns	20
S1235	Upper filter 35 µm for 2.5 ml columns - (thickness 1.5 mm)	20
S123532	Upper filter 35 µm for 2.5 ml columns - (thickness 3.2 mm)	20
S1290	Upper filter 90 µm for 2.5 ml columns	20
S1310	Upper filter 10 µm for 5.0 ml columns	20
S1335	Upper filter 35 µm for 5.0 ml columns - (thickness 1.5 mm)	20
S133532	Upper filter 35 µm for 5.0 ml columns - (thickness 3.2 mm)	20
S1390	Upper filter 90 µm for 5.0 ml columns	20
S1410	Upper filter 10 µm for 10 ml columns	20
S143532	Upper filter 35 µm for 10 ml columns - (thickness 3.2 mm)	20
S1490	Upper filter 90 µm for 10 ml columns	20
shipped at RT; store at RT		

### Column Sets

Order#	Product	Quantity
S10011	Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns (10 µm pore size)	3 each
S1001	Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns (35 µm pore size)	3 each
S10019	Set: Mobicol M1050; 2.5 ml; 5 ml and 10 ml columns (90 µm pore size)	3 each
S10031	Set upper filters 10 µm for S10011	1 set
S1003	Set upper filters 35 µm for S1001	1 set
S10039	Set upper filters 90 µm for S10019	1 set
shipped at RT; store at RT		

### Column Accessories

Order#	Product	Quantity
M3001	Tubing-to-Luer-lock adapter (material outside: metal; inside: glass and PTFE; with 1 m PTFE tubing) for Mobicols or laboratory columns	1
M3002	Luer adapter, (M+F) 10 each	20
M3003	Luer adapter for the top-cap (M)	20
M3004	Luer adapter for outlet (F)	20
M3006	Stopcocks, one-way, Luer-lock	4
M3007	Luer-lock adapter (M/M)	1
shipped at RT; store at RT		