SteMaxOne™ Media

Cell Culture Media for Human Mesenchymal Stem Cells

SteMaxOne™ is a chemically defined, protein and peptide-free culture media for human mesenchymal stem cells.


SteMaxOne™ medium can be used alone or as a basic nutrient mixture in combination with growth and maintenance factors selected by the user.



Available as liquid preparation ready for use or as liquid concentrate kits for higher volumes.

FEATURES:

- Free of hormones, proteins, peptides, and animal-derived components
- Free of complex additives such as albumines, hydrolysates, and yeast extract
- Exclusively made of small molecules listed in CAS- and EINECS registers
- Can be used as is or customized, unlike competing products



hMSC-TERT cells in SteMaxOne (Stem Cell 1) medium supplemented with 8 ng/mL bFGF and 2 mM L-glutamine.





www.bocascientific.com

Phone: (781) 686-1631 Toll free: (866) 227-6335

Email: info@bocascientific.com

