Product Information Sheet # 1042GE



SUMMARY shipped on blue ice; store at -20 °C For research use only. Not for use in diagnostic or therapeutic procedures.

Description

Polyclonal antisera were raised in rabbits after immunization with the conjugates: Citrullineglutaraldehyde carriers.

Contents: The antibody is supplied in a unit size of 50 µl in PBS with 50% glycerol.

Storage and Handling

When this product is stored undiluted at -20 °C or lower, it is stable for at least 2 years. Multiple freezing and thawing is only possible up to 5 times. The antibodies can be prediluted 10x in PBS containing 0.1% merthiolate or a mixture of PBS/glycerol (1:1, vol./vol.). At 4 °C these solutions are stable for 2 months.

Specificity

Using a conjugate Citrulline glutaraldehyde protein, antibody specificity was determined in an ELISA test by competition experiments with the following compounds:

| Compound | Cross-reactivity ratio ^(a) |
|------------------|---------------------------------------|
| Citrulline-G-BSA | 1 |
| Arginine-G-BSA | 1 : > 100,000 |
| Glutamine-G-BSA | 1 : > 100,000 |
| Ornithine-G-BSA | 1 : > 100,000 |

(a): Citrulline-G-BSA concentration / concentration of unconjugated or conjugated amino acids at half displacement ;

G = Glutaraldehyde, BSA = Bovine Serum Albumine

Recommended dilution

The anti-conjugated Citrulline antibodies can be diluted 1:2,000 to 1:5,000.

Product Information Sheet # 1042GE



Order Information, Shipping and Storage

| Order# | Product | Quantity |
|---------------|----------------------------|----------|
| 1042GE | Anti-Conjugated Citrulline | 50 µl |
| shipped on bl | ue ice; store at -20 °C | |

Contact and Support

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

| Customer Service – General inquiries & orders | Technical Service – Product information |
|--|--|
| phone: +49 (0)551 707 22 0 fax: +49 (0)551 707 22 22 e-mail: order@mobitec.com | phone: +49 (0)551 707 22 70 fax: +49 (0)551 707 22 77 e-mail: info@mobitec.com |
| MoBiTec in your area: Find your local distributor at | www.mobitec.com |