# **Anti-Conjugated NO-L-Glutamine**

Product Information Sheet # 1068GE



#### **SUMMARY**

shipped on blue ice; store at -20 °C

For research use only.

Not for use in diagnostic or therapeutic procedures.

#### **Description**

Polyclonal antisera were raised in rats after immunization with the conjugate: NO-L-Glutamine-glutaraldehyde-carriers.

**Contents:** The antibody is supplied in a unit size of 50 µl in PBS with 50% glycerol.

### Storage and Handling

When this product is stored undiluted at -20 °C or lower, it is stable for at least 2 years. Multiple freezing and thawing is only possible up to 5 times. The antibodies can be prediluted 10x in PBS containing 0.1% merthiolate or a mixture of PBS/glycerol (1:1, vol./vol.). At 4 °C these solutions are stable for 2 months.

### **Specificity**

Using a conjugate NO-L-Glutamine-glutaraldehyde protein, antibody specificity was determined in an ELISA test by competition experiments with the following compounds:

# Compound Cross-reactivity ratio<sup>(a)</sup>

(a): NO-L-Glutamine-G-BSA concentration / concentration of unconjugated or conjugated amino acids at half displacement;

G = Glutaraldehyde, BSA = Bovine Serum Albumine

#### **Recommended dilution**

The anti-conjugated NO-L-Glutamine antibodies can be diluted 1:2,000 to 1:5,000.

MoBiTec GmbH, Germany Phone: +49 551 70722 0 Fax: +49 551 70722 22 E-Mail: info@mobitec.com www.mobitec.com

# **Anti-Conjugated NO-L-Glutamine**

Product Information Sheet





## Order Information, Shipping and Storage

Order#	Product	Quantity
1068GE	Anti-Conjugated NO-L-Glutamine	50 µl
shipped on blue ice; store at -20 °C		

### **Contact and Support**

MoBiTec GmbH ● Lotzestrasse 22a ● D-37083 Goettingen ● Germany

Customer Service - General inquiries & orders Technical Service - Product information

phone: +49 (0)551 707 22 0 phone: +49 (0)551 707 22 70 fax: +49 (0)551 707 22 22 fax: +49 (0)551 707 22 77 e-mail: order@mobitec.com e-mail: info@mobitec.com

MoBiTec in your area: Find your local distributor at www.mobitec.com

MoBiTec GmbH, Germany Phone: +49 551 70722 0 Fax: +49 551 70722 22 E-Mail: info@mobitec.com www.mobitec.com