



Accessories for Affinity Chromatography Columns

The polypropylene Luer adapter (M) (M3003) for tubing connection to the inlet of the column allows a flow-through mode at low pressure, whereas the Luer adapter (F) (M3004) is suitable for a tubing connection to the outlet of the column. Furthermore, the Luer adapter (M/M) (M3007) facilitates the connection of the inlet to a Luer-lock system.

The Tubing-to-Luer-lock adapter (M3001) enables a safe connection from a pump to the inlet of a Mobicol or a 2.5 ml column.

The one-way stopcock, Luer-lock (F) (M3006) is a convenient tool for controlling flow-through at the in- or outlet.



M3003
Luer adapter (M) for Luer-lock cap



M3004
Luer adapter (F) for outlet



M3001 Tubing-to-Luer-lock adapter (connected to Mobicol "Classic" (M1002) with inserted large filter; not enclosed)



M3006
Stopcock, one-way, Luer-lock (F)



M3007
Luer-lock adapter (M/M)

Column Accessories

Order No.	Description	Amount
M3001	Tubing-to-Luer-lock adapter (material outside: metal; inside: glass and PTFE; with 1 m PTFE tubing) for Mobicols or laboratory columns	1
M3002	Luer adapter, (M+F) 10 each	20
M3003	Luer adapter for the top-cap (M)	20
M3004	Luer adapter for outlet (F)	20
M3007	Luer-lock adapter (M/M)	1
M3005	Bottom plugs for Mobicol M1002	50
M3006	Stopcocks, one-way, Luer-lock (F)	4
M3009	Luer-lock caps	50